# **TECHNOLOGY OFFER**

## **BLOOD MANAGER**


Blood Manager - Variable filling of hermetically sealed and preevacuated vials/containers via a pressure-adjusting adapter

#### BACKGROUND

Efficient use of blood and prevention of anemia are cornerstones of modern medicine. Patients with sepsis undergo repeated blood draws for blood culture testing. Blood is withdrawn into preevacuated and hermetically sealed vials. The vacuum in the vial is specified by the manufacturer and cannot be altered. Only 5-10 ml blood are necessary for blood culture testing. However, up to 30 ml blood are withdrawn without timely disconnection of the vial. Some patients require repeated determination of both aerobic and anaerobic samples (up to 3 pairs daily) leading to excess and preventable blood loss of 120 ml per day.

#### TECHNOLOGY

Our invention allows for specific withdrawal of a predefined sample volume (e.g. blood for blood culture testing) from a hermetically sealed and pre-evacuated vial/container. The adapter adjusts the pressure through a preset ancillary volume, and is suited for both aerobic and anaerobic samples



#### WISSENS/ TRANSFER/ OST/

www.wtz-ost.at

#### Tech ID: 816.19

#### **APPLICATIONS:** Blood cultures

DEVELOPMENT-STATUS: Prototype

**IPR:** AOF 586/2020 filed 05.05.2020

#### AVAIALABLE FOR:

Cooperation/licencing

#### **INVENTORS:**

David Baron (MedUniVienna) Christoph Eisenmenger-Sittner (TUVienna)

CONTACT: Christiane Galhaup PhD Medical University of Vienna

Technology Transfer Office Spitalgasse 23 A-1090 Vienna T: +43 40160 25206 christiane.galhaup@ meduniwien.ac.at www.meduniwien.ac.at



### ADVANTAGES

- Automatic stoppage of blood draw
- Defined sample volume
- Reduced blood loss